

# Tube Roller RM 2



## Description

Gentle tilting and rotating motion for uniform mixing

### Basic Equipment

Basic device with 6 removable plastic rolls and spill tray

### Loading Capacity

All standard tube sizes, flasks and universal vessels up to 80 mm in diameter..

## Details

|                          |                             |
|--------------------------|-----------------------------|
| <b>Motion:</b>           | rotating and tilting motion |
| <b>Shaking platform:</b> | 200 x 300 mm                |

|                                |                              |
|--------------------------------|------------------------------|
| <b>Max. load:</b>              | 2 kg                         |
| <b>Shaking speed:</b>          | 5 - 80 rpm, in steps of 5    |
| <b>Angle of tilt:</b>          | 3° fixed                     |
| <b>Stroke at end of rolls:</b> | 15 mm                        |
| <b>Dimensions of rolls:</b>    | 300 x Ø 32 mm                |
| <b>Number of rolls:</b>        | 6 pieces                     |
| <b>Runtime:</b>                | 0 - 99 h 59 min / continuous |
| <b>Order No.:</b>              | 6138 000                     |
| <b>Accessory:</b>              | Overhead Rack RM 2           |

## Technical specifications

|                                |                         |
|--------------------------------|-------------------------|
| <b>Electrical supply:</b>      | 100 - 240 V, 50 / 60 Hz |
| <b>Enclosure protection:</b>   | IP 21                   |
| <b>Heat emission:</b>          | max. 0,5 - 2 W          |
| <b>Ambient temperature:</b>    | 5°C to 45°C             |
| <b>Relative humidity:</b>      | ~ 85 % (non-condensing) |
| <b>Dimensions (w x d x h):</b> | 244 x 465 x 115 mm      |
| <b>Weight:</b>                 | 7,3 kg                  |