



FRONTIER™ 5706P

Multi Centrifuges



Reliable and Compact Multi-Function Low-Speed Centrifuge for General Purpose

Frontier 5706P Multi Centrifuge is a small benchtop multi-function centrifuge ideal for a multitude of separation applications. With versatile standard rotor options, the FC5706P seamlessly processes both round and conical bottom tubes, including common blood and urine tubes. Enhanced by user-friendly design, this centrifuge features an intuitive interface that streamlines operation and maximizes efficiency, all while conserving valuable bench space.

Unique Features:

- Designed for ease of operation, Frontier 5706P Multi Centrifuge features convenient control on an intuitive interface, which allows for responsive operation – even when wearing gloves.
- Enhanced adaptability to satisfy different customers' requirements in general laboratory applications. The FC5706P can hold up a variety of different tubes, starting from 50 ml to smaller volume 1,5 ml tubes. The Frontier 5706P offers 5 rotors options including a fixed angle rotor suitable for 18×1,5/2,0 ml tubes and a swing-out rotor that can hold 4×10 ml blood and urine tubes, boosting its versatility for various laboratory needs.
- Enclosing all its functionalities within a small footprint, the FC5706P saves valuable space on the lab bench.

FRONTIER™ 5706P *Multi Centrifuges*

Designed for a User-Friendly Experience

- Easy-to-use touch-wheel interface for parameter adjustments
- Convenient interface shows both set and actual parameters
- One-touch toggling between RPM/RCF and Acc/Dec
- Double-line LCD display with light blue back light
- Audible signal at the end of each run with a melody
- 10 acceleration/deceleration settings
- 99 memory sets to quickly store and recall commonly used parameters

Safety Features Enclosed Within Small Footprint

- Microprocessors control to help prevent accidents and protect both the user and the equipment
- Rotor imbalance sensor system stops the unit when it detects an imbalanced load
- Sturdy mechanical lid with locking system
- LED status to indicate the operational condition of the centrifuge like run completed (green), standby (yellow) or centrifugation (red)

Flexible Rotor Selections and Kits to Cover a Wide Range of General Laboratory and Diagnostic Applications

- Maintenance-free DC motor produces a maximum G-force of $4427 \times g$
- Speeds up to 6000 rpm and quick spin for short runs
- Additional fixed angle rotor ($18 \times 1,5/2,0$ ml microtubes) and swing out rotor (4×10 ml blood tubes) for enhanced flexibility




FRONTIER™ 5706P Multi Centrifuges

Specifications


Model		FC5706P
Speed Range		200 – 6 000 rpm; 10 rpm/set
Maximum RCF		4 427×g; 10 g/set
Maximum Capacity (Rotor)		Angle Rotor: 6×50 ml
Temperature range		Air cooled
Refrigerant		N/A
Running Time		Quick Spin, 10 sec to 99 hr 59 min, or continuous
Noise level (depending on the rotor)		≤ 56 ±2 dB(A)
Acc/Dec Levels		10/10
Allowable Density at Maximum Speed		1,2 g/ml
Allowable Kinetic Energy		2 427 Nm
Rated voltage		100-240 V/50-60 Hz
Voltage Fluctuation		± 10 %
Rated current		2,5 – 1,0 A
Power continues		70 W
Dimensions (W×D×H)		323×427×271 mm
Net Weight (without rotor)		14,2 kg
Shipping Dimensions (W×D×H)		370×540×365 mm
Shipping Weight (without rotor)		17,6 kg
Ambient conditions (EN/IEC 61010-1)	Environment	For indoor use only
	Altitude	Use up to an altitude of 2 000 m
	Ambient temperature	5 °C up to 35 °C
	Max. relative humidity	Max. relative humidity 80 % for temperatures up to 31 °C Decreasing linearly to 50 % relative humidity up to 35 °C.
EMC		EN/IEC 61326-1 Class B emissions, Basic immunity FCC Class B emissions
Not suitable for use in hazardous environments.		

Accessories

	Order No.	30130877				
	Description	Rotor Angle 12×15 ml RB/FA				
	Included in Kit	31012113				
	Max. Speed	6 000 rpm				
	Max. RCF	4 427×g				
Accessories						
Order No.	Description	Vol./Tube	Diameter	Tube Type	No.of Tubes/ Rotor	
Without Adapter	Rotor Angle 12×15 ml RB/FA	15 ml	17 mm	Round Bottom/Conical Bottom	12	
30211180	Adapter 1×6-10 ml D16 mm Vac/Sar, 2 pk	6-10 ml	16 mm	Vacutainer®/Sarsted	12	
30130889	Adapter 1×4-7 ml D13 mm Vac/Sar, 2 pk	4-7 ml	13 mm	Round Bottom/Vacutainer®/Sarstedt	12	
30130890	Adapter 1×1,8-5 ml D13 mm Vac/Sar, 2 pk	1,8-5 ml	13 mm	Round Bottom/Vacutainer®/Sarstedt	12	
30130886	Adapter 1×1,5/2,0 ml D11mm, 6 pk	1,5/2,0 ml	11 mm	Microtube	12	

FRONTIER™ 5706P Multi Centrifuges

Accessories

	Order No.	30130878
	Description	Rotor Angle 6 × 50 ml RB/FA
	Included in Kit	31012115
	Max. Speed	6 000 rpm
	Max. RCF	4 427 × g

Accessories

Order No.	Description	Vol./Tube	Diameter	Tube Type	No.of Tubes/ Rotor
Without Adapter	Rotor Angle 6 × 50 ml RB/FA	50 ml	29 mm	Round Bottom/Conical Bottom	6
30130891	Adapter 1 × 30 ml D26 mm, 2 pk	30 ml	26 mm	Round Bottom	6
30130892	Adapter 1 × 16 ml D18 mm, 2 pk	16 ml	18 mm	Round Bottom	6
30130893	Adapter 1 × 15 ml D17 mm, 2 pk	15 ml	17 mm	Round Bottom	6
30130894	Adapter 1 × 15 ml D17 mm RB/FA, 2 pk	15 ml	17 mm	Round Bottom/Conical Bottom	6

	Order No.	30130880
	Description	Rotor Swing out 6 × 5 ml RB
	Max. Speed	4 000 rpm
	Max. RCF	1 878 × g


Accessories

Order No.	Description	Vol./Tube	Diameter	Tube Type	No.of Tubes/ Rotor
30130895	Buckets For 30130880, 2 pk*	5 ml	13 mm	Round Bottom/Vacutainer®/Sarstedt	6

	Order No.	31013397
	Description	Rotor Swing Out 4 × 10 ml
	Included in Kit	31012114
	Max. Speed	4 200 rpm
	Max. RCF	2 406 × g

Accessories

Order No.	Description	Vol./Tube	Diameter	Tube Type	No.of Tubes/ Rotor
31021231	Rack 1 × 1,8-10 ml D16 mm Vac/Sar, 4 pk**	1,8-10 ml	16 mm	Round Bottom/ Vacutainer®/ Sarstedt	4

	Order No.	83041499
	Description	Rotor Angle 18 × 1,5/2,0 ml
	Max. Speed	6 000 rpm
	Max. RCF	2 938 × g

Accessories

Order No.	Description	Vol./Tube	Diameter	Tube type	No.of Tubes/ Rotor
Without Adapter	Rotor Angle 18 × 1,5/2,0 ml	1,5/2,0 ml	11 mm	Microtube	18
30130884	Adapter 1 × 0,2/0,4 ml D6 mm, 6 pk	0,2/0,4 ml	6 mm	Microtube	18
30130885	Adapter 1 × 0,5 ml D8 mm, 6 pk	0,5 ml	8 mm	Microtube	18

*Buckets included

**Racks included

Shortcut list:

Vac – Vacutainer® Sar – Sarstedt FA – Conical bottom / Falcon® FB – Flat bottom RB – Round bottom

FRONTIER™ 5706P *Multi Centrifuges*


Other Standard Features and Equipment

Microprocessor control, Imbalance switch-off, Audible signal at the end of each run with melody options, Emergency release, 99 memory sets, Ten acceleration and ten deceleration rates, Pre-selection of running time from 10 sec to 99 hr 59 min or continuous, Quick spin for short runs.

Compliance

See Declaration of Conformity available CE Declarations | OHAUS

Outline Dimensions



Click or scan
for the web



80777081_A 20250418 © Copyright OHAUS Corporation

OHAUS Europe GmbH

e-mail: ssc@ohaus.com

Tel: 0041 22 567 53 19

e-mail: tsc@ohaus.com

Tel: 0041 22 567 53 20

www.ohaus.com

*The management system
governing the manufacture
of this product is
ISO 9001:2015 certified.*

