


A 100

Upper basket

For the optimum loading of injector modules or inserts.



- Upper basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Loading height with module with plastic supports 190 mm
- Clearances: width 489, length 440 mm

EAN: 4002516841296 / Material number: 12642640 / Old Material Number: 69610002D

Product category	
Washer-disinfectors, laboratory	•
Cleaning machines for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
Product category	
Basket with module connection	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	141
External dimensions, net width in mm	525
External dimensions, net depth in mm	528
External dimensions, gross height in mm	170
External dimensions, gross width in mm	555
External dimensions, gross depth in mm	565
Net weight in kg	2.92
Gross weight in kg	3.56