


A 101

Upper basket

Height adjustable, for the optimum loading of inserts.



- Vertical clearance 160 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm


EAN: 4002516828143 / Material number: 12645590 / Old Material Number: 69610101D

Product category	
Fresh-water dishwashers	•
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Cleaning machines for industrial parts cleaning	•
Washer-disinfectors, fire brigades and emergency services	•
Machine affiliation	
PG 8063	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Upper basket	•
Carrier/basket carrier	•
Upper basket	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of different utensils	•
Reprocessing of breathing air regulators	•
Dimensions and weight	
External dimensions, net height in mm	206
External dimensions, net width in mm	528
External dimensions, net depth in mm	527
External dimensions, gross height in mm	220
External dimensions, gross width in mm	540

A 101

Upper basket

Height adjustable, for the optimum loading of inserts.



EAN: 4002516828143 / Material number: 12645590 / Old Material Number: 69610101D

External dimensions, gross depth in mm	550
Net weight in kg	3.03
Gross weight in kg	3.68