


A 102

Upper basket

Height adjustable, for the optimum loading of inserts.



- Vertical clearance 205 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm


EAN: 4002516855262 / Material number: 12645600 / Old Material Number: 69610201D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Cleaning machines for industrial parts cleaning	•
<b>Machine affiliation</b>	
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Carrier/basket carrier	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
<b>Dimensions and weight</b>	
External dimensions, net height in mm	206
External dimensions, net width in mm	528
External dimensions, net depth in mm	527
External dimensions, gross height in mm	230
External dimensions, gross width in mm	540

A 102

Upper basket

Height adjustable, for the optimum loading of inserts.



EAN: 4002516855262 / Material number: 12645600 / Old Material Number: 69610201D

External dimensions, gross depth in mm	545
Net weight in kg	3.04
Gross weight in kg	3.78