


A 14/1
Lid
for securing delicate load items.




- Lid in stainless steel
- 8 x 8 mm perforations, 3 mm slat
- Covers half of AK 12/1 insert

EAN: 4002516340126 / Material number: 11131560 / Old Material Number: 69746003D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
PLW 8505	•
Product category	
Cover	•
Lid	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•

A 14/1
Lid
for securing delicate load items.



EAN: 4002516340126 / Material number: 11131560 / Old Material Number: 69746003D

Dimensions and weight	
External dimensions, net height in mm	19
External dimensions, net width in mm	207
External dimensions, net depth in mm	224
External dimensions, gross height in mm	65
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	250
Net weight in kg	0.20
Gross weight in kg	0.35