

A 151

Lower basket

For holding standard DIN mesh trays as well as various inserts.



- Vertical clearance dependent on upper basket selected
- Loading width: 490 mm, loading depth: 498 mm

EAN: 4002516853992 / Material number: 12666340 / Old Material Number: 69615101D

Product category	
Fresh-water dishwashers	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Cleaning machines for industrial parts cleaning	•
Washer-disinfectors, fire brigades and emergency services	•
Machine affiliation	
PG 8058	•
PG 8063	•
PLW 8636 [LAB / LAB MON]	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Carrier/basket carrier	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of mesh trays	•

A 151

Lower basket

For holding standard DIN mesh trays as well as various inserts.



EAN: 4002516853992 / Material number: 12666340 / Old Material Number: 69615101D

Reprocessing of different utensils	•
Reprocessing of breathing air regulators	•
Reprocessing of breathing masks	•
Reprocessing of compressed air unit carrier frames	•
<b>Capacity</b>	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
<b>Dimensions and weight</b>	
External dimensions, net height in mm	88
External dimensions, net width in mm	529
External dimensions, net depth in mm	522
External dimensions, gross height in mm	120
External dimensions, gross width in mm	545
External dimensions, gross depth in mm	535
Net weight in kg	1.84
Gross weight in kg	2.40