


A 200
Mobile unit
for the optimum loading of 38 pipettes.



- for 10 100 ml pipettes in row 1, distance 20 mm
- for 14 25 ml pipettes in a 2nd row, distance 26 mm
- for 14 10 ml pipettes in a 3rd row, distance 26 mm

EAN: 4002516854012 / Material number: 12667060 / Old Material Number: 69620001D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
Product category	
Mobile unit for pipettes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Pipettes [number]	38
Dimensions and weight	
External dimensions, net height in mm	397
External dimensions, net width in mm	529
External dimensions, net depth in mm	546
External dimensions, gross height in mm	520
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	550
Net weight in kg	4.95
Gross weight in kg	6.38