


A 301/5

Injector module

for optimum loading of laboratory glassware – can be equipped with 18 nozzles.

EAN: 4002516188971 / Material number: 11116950 / Old Material Number: 69630105D



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Equipped with 9 A 842 nozzles and 9 A 843 nozzles
- Part of the EasyLoad system
- Nozzle spacing 78 mm across the width, 78 mm across the depth

Product category	
Washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	18
Dimensions and weight	
External dimensions, net height in mm	203
External dimensions, net width in mm	232
External dimensions, net depth in mm	480
External dimensions, gross height in mm	200
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	515
Net weight in kg	2.16
Gross weight in kg	3.17
Companion pack	
A 842, injector nozzle Ø 4 mm, length 90 mm [number]	9
A 843, injector nozzle Ø 4 mm, length 185 mm [number]	9