

A 304 Module

For the optimum loading of up to 98 vials.



- Suitable for 98 centrifuge tubes, phials, test tubes, for example
- The loading heights differ depending on the positioning
- 1 unit can be inserted into each of the baskets A 100, A 150 and APLW 106

EAN: 4002516854982 / Material number: 12669230 /
Old Material Number: 69630401D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of test tubes	•
Capacity	
Phials [number]	98
Dimensions and weight	
External dimensions, net height in mm	130
External dimensions, net width in mm	222
External dimensions, net depth in mm	471
External dimensions, gross height in mm	150
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	490
Net weight in kg	2.30
Gross weight in kg	3.15