

A 306/1 Module

for the optimum loading of large-volume laboratory glassware.



- Suitable for 4 1–2 litre measuring cylinders, high shape
- Suitable for 2 1–2 litre measuring cyl., low shape & 2 high shape
- Plastic-coated contact surfaces
- Measuring cylinders even removable when module is in machine
- Maximum length of the measuring cylinders in the A 150 lower basket: 500 mm
- Up to 2 modules can be inserted in the baskets A 150 and APLW 106

EAN: 4002516854456 / Material number: 12683100 /
Old Material Number: 69630602D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Insert/module for measuring cylinder	•
Module for large-volume laboratory glassware	•
Product properties	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Colour	Stainless Steel, Black, White
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Capacity	
Measuring cylinders [number]	4
Dimensions and weight	
External dimensions, net height in mm	235
External dimensions, net width in mm	418
External dimensions, net depth in mm	471
External dimensions, gross height in mm	250
External dimensions, gross width in mm	390
External dimensions, gross depth in mm	485
Net weight in kg	2.87
Gross weight in kg	3.41