

A 312
Module
For reprocessing viscometers.



- Suitable for up to 8.0 viscometers
- Equipped with 8 x A 867 special nozzles
- Nozzle spacing 118 mm across the width, 118 mm across the depth

EAN: 4002516087397 / Material number: 11054470 / Old Material Number: 69631200D

Product category	
Washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Module for viscometers	•
Product properties	
Material	Stainless steel, Plastic, Silicone
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Viscometer [quantity]	8.0
Dimensions and weight	
External dimensions, net height in mm	66
External dimensions, net width in mm	133
External dimensions, net depth in mm	479
External dimensions, gross height in mm	250
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	505
Net weight in kg	1.58
Gross weight in kg	2.33