

A 313

Module

For reprocessing large-volume laboratory glassware.



- Laboratory flasks, up to 10 l
- Round flasks, up to 6 l
- Erlenmeyer flasks, up to 5 l
- Measuring flasks, up to 5 l
- Can be equipped with nozzles & holders based on individual needs
- Equipped with 5 injector positions for 2 or 3 load items

EAN: 4002516087854 / Material number: 11055970 / Old Material Number: 69631300D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
<b>Machine affiliation</b>	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
<b>Product category</b>	
Insert/module for narrow-necked glassware	•
Module for large-volume laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	221
External dimensions, net width in mm	234
External dimensions, net depth in mm	495
External dimensions, gross height in mm	260
External dimensions, gross width in mm	270
External dimensions, gross depth in mm	505
Net weight in kg	2.82
Gross weight in kg	3.29