

A 846

Injector nozzle

with support, for cleaning laboratory flasks and round flasks.



- For laboratory bottles, max. 10 l, and round flasks, max. 6 l
- Nozzle length: 127 mm, nozzle length with holder: 158 mm
- Nozzle Ø: 10 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516855088 / Material number: 12665100 / Old Material Number: 69684601D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
<b>Machine affiliation</b>	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Rubber, Plastic
Colour	Stainless Steel, Black, Grey
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	158
External dimensions, net width in mm	110
External dimensions, net depth in mm	110
External dimensions, gross height in mm	160
External dimensions, gross width in mm	145
External dimensions, gross depth in mm	265
Net weight in kg	0.23
Gross weight in kg	0.38