


A 847

Injector nozzle
with support, for cleaning Erlenmeyer flasks.



- For Erlenmeyer flasks, max. 5 l
- Nozzle length: 127 mm
- Nozzle Ø: 10 mm, Ø of support: 142 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516855101 / Material number: 12665110 / Old Material Number: 69684701D

Product category	
Washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	170
External dimensions, net width in mm	249
External dimensions, net depth in mm	35
External dimensions, gross height in mm	70
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	260
Net weight in kg	0.13
Gross weight in kg	0.22