

A 848

Injector nozzle  
with support, for cleaning measuring cylinders.



- For measuring flasks, max. 5 l
- Nozzle length: 350 mm
- Nozzle Ø: 10 mm, Ø of support: 142 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516854234 / Material number: 12665120 / Old Material Number: 69684801D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
<b>Machine affiliation</b>	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	350
External dimensions, net width in mm	20
External dimensions, net depth in mm	20
External dimensions, gross height in mm	70
External dimensions, gross width in mm	185
External dimensions, gross depth in mm	360
Net weight in kg	0.20
Gross weight in kg	0.32