


A 862
Holder rack
for A 302/2 and A 302/3 injector modules.



- Arrangement: 4 x 8 fields
- For use on A 302/2 and A 302/3 injector modules
- For the optimum positioning of laboratory glassware
- Plastic-coated

EAN: 4002516854203 / Material number: 12665480 /
Old Material Number: 69686201D

Product category	
Washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Other inserts/other modules	•
Product properties	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Colour	Stainless Steel, Black
Number of compartments	32
Compartment width in mm	56
Compartment length in mm	52
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	30
External dimensions, net width in mm	250
External dimensions, net depth in mm	350
External dimensions, gross height in mm	100
External dimensions, gross width in mm	255
External dimensions, gross depth in mm	465
Net weight in kg	0.55
Gross weight in kg	0.80