

ADB
Cover plate
For built-under versions.



- To protect the worktop from moisture
- Optional when installed as a built-under appliance

EAN: 4002515379493 / Material number: 09715590 / Old Material Number: 67216008D

Product category	
Fresh-water dishwashers	•
Machine affiliation	
PG 8055 U	•
PG 8056 U	•
PG 8057 TD U	•
PG 8058 U	•
PG 8059 U	•
PG 8061 U	•
PG 8096 U	•
PG 8099 U	•
PFD 400 U	•
PFD 401 U	•
PFD 402 U	•
PFD 404 U	•
PFD 405 U	•
PFD 407 U	•
PFD 408 U	•
PG 8582	•
PG 8592	•
Product category	
Accessories for built-under machines	•
Installation materials	•
Product properties	
Material	Stainless steel, Galvanised sheet steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	25
External dimensions, net width in mm	596
External dimensions, net depth in mm	542
External dimensions, gross height in mm	70
External dimensions, gross width in mm	630
External dimensions, gross depth in mm	650
Net weight in kg	2.69
Gross weight in kg	3.75