

DTD 2  
Drip tray  
for the secure placement of ProCare media canisters.



- Space for 2 ProCare media canisters
- Secure place for 5 and 10 l canisters in the ProCare series
- Optimum residual emptying of the canister due to angle
- Suitable for use in an adjacent unit

EAN: 4002515933695 / Material number: 10695460 / Old Material Number: 61910001D

<b>Product category</b>	
ProCare	•
<b>Applications</b>	
Secure placement of ProCare media canisters	•
Optimum emptying of ProCare media canisters	•
<b>Machine affiliation</b>	
PFD 100 (SmartBiz)	•
PFD 101	•
PFD 101 i	•
PFD 101 U	•
PFD 102 i	•
PFD 103 SCi XXL	•
PFD 104 SCVi XXL	•
PG 8063	•
PTD 901	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
<b>Product category</b>	
Other accessories	•
<b>Product properties</b>	
Material	Plastic
Colour	Grey
Option to accommodate 2 ProCare media canisters 5 or 10 litres in size	•
Inclined floor for optimum residual emptying of ProCare media canisters	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	62
External dimensions, net width in mm	262
External dimensions, net depth in mm	420
External dimensions, gross height in mm	74.00
External dimensions, gross width in mm	327.00
External dimensions, gross depth in mm	463.00

DTD 2  
Drip tray  
for the secure placement of ProCare media canisters.



EAN: 4002515933695 / Material number: 10695460 / Old Material Number: 61910001D

Gross weight in kg	0.72
<b>Capacity</b>	
Canisters, maximum size 190 x 230 mm [number]	2
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	74
External dimensions, gross width in mm	327
External dimensions, gross depth in mm	463
Gross weight in kg	1
<b>Dimensions and weight</b>	
External dimensions, net height in mm	62
External dimensions, net width in mm	420
External dimensions, net depth in mm	262
External dimensions, gross height in mm	74
External dimensions, gross width in mm	327
External dimensions, gross depth in mm	463
Net weight in kg	0.42
Gross weight in kg	0.72