


E 103/1

Insert

For the optimum loading of test tubes up to 12 x 75 mm.



- Suitable for about 200 Test tubes (up to 12 x 75 mm)
- Divided into 6 surfaces
- The mesh width is 8 x 8 mm
- Including lid A 13

EAN: 4002516854104 / Material number: 12667340 / Old Material Number: 69510303D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	200
Dimensions and weight	
External dimensions, net height in mm	122
External dimensions, net width in mm	200
External dimensions, net depth in mm	220
External dimensions, gross height in mm	140
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	250
Net weight in kg	0.56
Gross weight in kg	0.73
Companion pack	
A 13 lid [number]	1

E 103/1

Insert

For the optimum loading of test tubes up to 12 x 75 mm.