


E 105/1

Insert

For the optimum loading of test tubes up to 12 x 165 mm.



- Suitable for about 200 Test tubes (up to 12 x 165 mm)
- Capacity including lid A 13
- The mesh width is 9 x 9 mm

EAN: 4002516853978 / Material number: 12667380 / Old Material Number: 69510503D

|  |                 |
|--|-----------------|
| <b>Product category</b>  |                 |
| Washer-disinfectors, laboratory                                | •               |
| Large-capacity washer-disinfectors, laboratory                 | •               |
| SlimLine laboratory washers                                    | •               |
| Cleaning machines for industrial parts cleaning                | •               |
| Large-chamber laboratory washers for industrial parts cleaning | •               |
| SlimLine laboratory washers for industrial parts cleaning      | •               |
| <b>Machine affiliation</b>                                     |                 |
| G 7825   | •               |
| PG 8527  | •               |
| PG 8535  | •               |
| PG 8536  | •               |
| PLW 8636 [LAB / LAB MON]                                       | •               |
| PG 8504  | •               |
| PG 8583  | •               |
| PG 8583 CD   | •               |
| PG 8593  | •               |
| PLW 8615   | •               |
| PLW 8616   | •               |
| PLW 8617   | •               |
| PLW 8683   | •               |
| PLW 8683 CD  | •               |
| PLW 8693   | •               |
| PLW 6011   | •               |
| PLW 6111   | •               |
| PLW 7111   | •               |
| PLW 8505   | •               |
| <b>Product category</b>  |                 |
| Insert/module for test tubes                                   | •               |
| <b>Product properties</b>                                      |                 |
| Material   | Stainless steel |
| Colour   | Stainless Steel |
| <b>Area of application</b>                                     |                 |
| Reprocessing of test tubes                                     | •               |
| Reprocessing of laboratory glassware                           | •               |
| <b>Capacity</b>  |                 |
| Test tubes [number]  | 200             |
| <b>Dimensions and weight</b>                                   |                 |
| External dimensions, net height in mm                          | 212             |
| External dimensions, net width in mm                           | 200             |
| External dimensions, net depth in mm                           | 210             |
| External dimensions, gross height in mm                        | 235             |
| External dimensions, gross width in mm                         | 240             |
| External dimensions, gross depth in mm                         | 250             |
| Net weight in kg   | 0.88            |
| Gross weight in kg   | 1.15            |
| <b>Companion pack</b>  |                 |
| A 13 lid [number]  | 1               |

E 105/1

Insert

For the optimum loading of test tubes up to 12 x 165 mm.