

E 139/1

Insert

For the optimum loading of test tubes up to 12 x 200 mm.



- Suitable for test tubes up to 12 x 200 mm
- Divided into 6 surfaces
- Mesh width 9 x 9 mm
- For use in lower basket
- Including lid A 13

EAN: 4002516854036 / Material number: 12667390 / Old Material Number: 69513903D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	200
<b>Dimensions and weight</b>	
External dimensions, net height in mm	243
External dimensions, net width in mm	200
External dimensions, net depth in mm	230
External dimensions, gross height in mm	245
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	265
Net weight in kg	1.00
Gross weight in kg	1.21
<b>Companion pack</b>	
A 13 lid [number]	1

E 139/1

Insert

For the optimum loading of test tubes up to 12 x 200 mm.