

ProCare Lab 30 C - 5 l

Neutralising agent, acidic, 5 l

for material-saving machine cleaning of laboratory glassware and utensils.



- Free from surfactants and phosphates
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics and plastics
- To be used as neutralising agent or acid pre-cleaner
- Removes residues that are soluble in acids from sensitive materials

EAN: 4002515660744 / Material number: 10266670 /
Old Material Number: 62927101EU1

Product category	
Washer-disinfectors, laboratory	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing	
pH value at 20 °C	1,0
Dispensing recommendation	
Dispensing recommendation per wash cycle	1-4 ml/l (0.1-0.4%)
Dispensing recommendation (per cycle)	1-4 ml/l (0.1-0.4%)
Large opening for suction lance	•
Suitable for all levels of water hardness	•
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Industrial laboratories	•
Phosphate industry	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Acidic neutralising agent
Aggregate form	Liquid
Ingredients	Organic acids
Minimum storage temperature in °C	-3
Maximum storage temperature in °C	30
Dimensions and weight	
External dimensions, net height in mm	217
External dimensions, net width in mm	142
External dimensions, net depth in mm	248
External dimensions, gross height in mm	250.00
External dimensions, gross width in mm	190.00
External dimensions, gross depth in mm	150.00
Gross weight in kg	6.24
Fill volume in ml	5000
Product properties	
Medium type	Acidic neutralising agent
Aggregate form	Liquid
Composition	Organic acids
Minimum storage temperature in °C	-3
Maximum storage temperature in °C	30
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Suitable for sensitive materials	•

ProCare Lab 30 C - 5 l

Neutralising agent, acidic, 5 l

for material-saving machine cleaning of laboratory glassware and utensils.



EAN: 4002515660744 / Material number: 10266670 /

Old Material Number: 62927101EU1

Suitable for neutralisation of non-anodised aluminium	•
Not suitable for anodised aluminium	•
Not suitable for light and non-ferrous metals	•
Not suitable for chrome-plated parts	•
Not suitable for nickel-plated parts	•
Not suitable for acidic pre-cleaning of non-anodised aluminium	•
Does not contain phosphates	•
Does not contain surfactants	•
Very gentle on materials	•
Acidic pre-cleaning	•
Neutralisation of carried-over alkali residues	•
Dimensions and weight	
External dimensions, net height in mm	217
External dimensions, net width in mm	142
External dimensions, net depth in mm	248
External dimensions, gross height in mm	250
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	150
Gross weight in kg	6
Fill volume in ml	5000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•