


### Multi-Tube Vortex Mixer



**Features:**

- \* Mixing mode changes quickly.
- \* LCD Integrated display.
- \* High speed accuracy.
- \* Suitable for a variety of test tube racks.
- \* The appearance and operation interface are simple and beautiful.

**Picture:**



**Technical Parameters:**

Model	BK-VX2800
Speed Range	300-2800rpm
Speed Regulation Accuracy	±1%
Amplitude	3mm
Time Setting Range	1s-99h59min
Pulse Interval Timing Range	1s-99s
Pulse Operation Timing Range	1s-999s
Max Loading	5.5kg
Consumption	70W
Power Supply	AC100-240V, 50/60Hz
Top Cover Size	310*210mm
External Size (L*W*H)	420*305*465mm
Net Weight	20kg
Package Size (W*D*H)	492*378*520mm
Gross Weight	23kg

**Accessories:**

Model	Specification	Compatible Tube Type	Hole NO.	Size
<b>Standard</b>				
A1	Replaceable tray mat	/	/	305*204*50mm
A2	Φ17mm plastic test tube rack	96*10ml virus sampling tube can be placed	96	294*200*50mm
A3	Φ26mm plastic test tube rack	Can put 20 mixed 1 virus sampling tube rack	48	280*214*53mm
A4	Φ16mm plastic test tube rack	Can put 10ml/15ml centrifugal tube	50	245*132*45mm
A5	Φ29mm plastic test tube rack	50ml centrifuge tube	15	245*132*45mm
<b>Optional</b>				
A6	Φ10mm(Φ12mm) sponge test tube rack	/	50	245*132*45mm
A7	Φ13mm sponge test tube rack	/	50	245*132*45mm
A8	Φ13mm sponge test tube rack	/	15	245*132*45mm
A9	Φ13mm sponge test tube rack	/	10	245*132*45mm